Verify*Now*[®] IIb/IIIa

PLATELET REACTIVITY TEST

Accu:metrics[®]

INSTRUCTIONS FOR USE

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Technical Support
(US) +1 858 643 1600
support@accumetrics.com

INTENDED USE

The VerifyNow IIb/IIIa Test is a semi-quantitative, whole blood platelet function test used to measure glycoprotein (GP) IIb/IIIa receptor blockade in patients treated with abciximab or eptifibatide. VerifyNow IIb/IIIa Test results should be interpreted in conjunction with other clinical and laboratory data available to the clinician.

For in vitro diagnostic use.

SUMMARY AND EXPLANATION

The VerifyNow System is a turbidimetric based optical detection system, which measures platelet-induced aggregation. The system consists of an instrument, a disposable test device and quality control materials. See **Figure 1** for a representation of the test device. Quality control measures include an instrument based electronic quality control (EQC), two levels of wet quality controls (WQC), internal quality controls, and shipping controls. The instrument controls all assay test sequencing, temperature, reagent sample mixing and performs self-diagnostics. The degree of platelet function is determined and the result is displayed.


Figure 1: Test Device

Each single use test device contains a lyophilized preparation of human fibrinogen coated beads, thrombin receptor activating peptide (iso TRAP), and buffer. The patient sample is anticoagulated whole blood, which is automatically dispensed from the blood collection tube into the test device by the instrument, with no blood handling required by the user.

PRINCIPLE OF THE TEST

The VerifyNow IIb/IIIa Test is designed to assess platelet function based upon the ability of activated platelets to bind fibrinogen. Fibrinogen-coated microparticles aggregate in whole blood in proportion to the number of unblocked platelet GPIIb/IIIa receptors. The rate of microbead aggregation is more rapid and reproducible if platelets are activated; therefore the reagent iso-TRAP is incorporated into the test to induce platelet activation without fibrin formation. Light transmittance increases as activated platelets bind and aggregate fibrinogen-coated beads. The instrument measures this change in optical signal and reports results in Platelet Aggregation Units (PAU).

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use.

The VerifyNow Instrument and its components should only be used as directed in the User Manual.

Do not use the VerifyNow IIb/IIIa Test device or WQC materials beyond the expiration date.

All patient samples should be handled as if capable of transmitting disease.

The reagents are manufactured with a material purified from human plasma that was found negative for all communicable diseases tested, including HIV-1, HIV-2, Hepatitis B surface antigen (HBsAg) and HCV. Handle test devices as biohazardous material and dispose of in an appropriate manner.

Do not attempt to remove the blood collection (or WQC) tube from the test device after it has been inserted onto the needle in the test device.

REAGENT STORAGE AND HANDLING

Test Device Kit Indicator: Each VerifyNow Test kit has a temperature indicator on the outside of the packaging. The user is instructed to inspect the indicator upon receipt of the kit. If the indicator has changed color, the kit has been exposed to elevated temperature, and a Wet Quality Control (WQC) Level 2 must be run to ensure that the reagents are performing properly. The Temperature Indicator detects errors due to adverse environmental conditions.

For extended storage, maintain test devices at 2° to 8°C. Alternatively, these test devices may be stored for up to 8 weeks at room temperature (18° to 25°C) or (64° to 77°F), but not beyond the expiration date printed on the pouch. When removed from refrigeration, write the date to discard on the pouches and/or kit box.

Allow test devices to reach room temperature prior to use.

Test devices should remain sealed in the foil pouches until ready for use to prevent damage by humidity.

Test devices are stable until the expiration date printed on the outer box.

SAMPLE COLLECTION PROCEDURE

Instructions for Blood Collection Tube Selection for All Samples

Samples must be collected into the appropriate sample tube and anticoagulant according to the table below. Samples collected from patients who are receiving eptifibatide (Integrilin) must only be collected into tubes containing lithium heparin.

Greiner 454322 (2 mL citrate tube)	Greiner 454323 (high altitude 2mL citrate tube)	Greiner 454244 (3 mL lithium heparin tube)
---------------------------------------	--	---

Table 1: Blood Collection Tubes

Instructions for Sample Collection from Indwelling Venous Catheters:

Note: Do not use an Intravenous Catheter

1. Whole blood samples that are obtained from an indwelling catheter should be collected after sufficient discard (approximately 5 mL) has been drawn to clear the line. Ensure indwelling catheter is free of clots. **Transfer blood to the blood collection tube immediately after collection.**

2. Gently invert the sample tube **5 times** to ensure complete mixing of the contents.

3. Test the sample within 15 minutes of collection as instructed in the Test Procedure.

4. Blood drawn during the insertion of an IV catheter should not be used.

Instructions for Peripheral Samples:

1. Whole blood may be collected from venous sites using a 21 gauge or larger (e.g. 18-20 gauge) needle in an appropriate blood collection tube. Blood samples should be obtained from an extremity free of peripheral venous infusions.

2. Collect a discard tube first (approximately 2 mL). The discard tube must not contain EDTA.

3. Gently invert the sample tube at least 5 times to ensure complete mixing of the contents.

4. Test the sample within 15 minutes of collection as instructed in the Test Procedure.

5. Blood drawn during the insertion of an IV catheter should not be used.

SAMPLE COLLECTION PRECAUTIONS

- Improper blood collection techniques may cause an Attention message. Refer to the VerifyNow System User Manual for a more detailed explanation of this message.
- Use only 21 gauge or larger bore needles for blood collection or transfer.
- Blood samples should be kept upright prior to testing and avoid prolonged contact with the rubber stopper on the blood collection tube.
- Avoid use of a rocker or pneumatic tube transport system.
- Collection of the blood sample must be performed with care to avoid hemolysis or contamination by tissue factors. Samples with evidence of clotting should not be used.
- The first collection tube must be discarded (approximately 2 mL).
- If drawing blood for a CBC at the same time as sample collection for VerifyNow IIb/IIIa Test, draw the VerifyNow sample first. Fill the CBC tube after the VerifyNow sample tube.
- Fresh whole blood samples must be used within 15 minutes of collection.
- Always ensure blood collection tubes are filled to the indicated fill volumes. At altitudes greater than 2500 feet above sea level, blood collection tubes may not fill to the specified volume, which results in an incorrect ration of blood to anticoagulant. Users at these altitudes should refer to their facility's blood collection protocols or blood collection tube manufacturer's recommendations for instructions to properly fill blood collection tubes.
- Do not freeze or refrigerate blood samples.
- All patient samples should be handled as if capable of transmitting disease.
- Universal precautions should be followed.

TEST PROCEDURE

1. Select the type of sample that will be tested: VerifyNow WQC or a patient blood sample.

2. Open the foil pouch and remove the test device. Test devices should only be handled by finger grip (**See Figure 1**).

3. Remove the needle's protective sheath by pulling directly up on the sheath. Do not twist the sheath as this may remove the needle.

4. At the instrument prompt, insert the test device into the instrument.

5. At the instrument prompt, invert the sample tube 5 times, and insert onto the needle in the test device. Close the test port cover.

6. The instrument will run the test and display the result.

The sample is under pressure. Do not remove the sample tube from the test device.

7. Remove the test device by grasping the test device finger grip and pulling straight up. Do not remove the tube from the test device. Dispose of the entire test device/sample tube in appropriate biohazard waste container.

Refer to the VerifyNow System User Manual for complete operating instructions.

MATERIALS PROVIDED

10 VerifyNow IIb/IIIa Test devices (Catalog # 85310) individually sealed in foil pouches. Each test device contains lyophilized fibrinogen-coated beads, iso-TRAP, bovine serum albumin, and buffer.

MATERIALS REQUIRED BUT NOT PROVIDED

- VerifyNow System (Catalog #86000)
- VerifyNow WQC Kit (Catalog #85047)
- Blood collection tubes listed in **Table 1**.

Note: Sodium citrate should not be used as an anticoagulant when measuring samples from patients treated with eptifibatide.

QUALITY CONTROL

Instrument electronic controls

- To ensure instrument performance, the manufacturer recommends that an Electronic Quality Control (EQC) be run once each day of patient testing. This reusable device verifies instrument optics, pneumatics and mixing. Refer to the Quality Control section of the VerifyNow System User Manual for instructions on running EQC.
- The instrument automatically verifies sample filling, correct fluid transfer, and mixing. The instrument also monitors the electronic and mechanical components.

Test device controls

The VerifyNow IIb/IIIa Test contains the following internal control:

- **Negative control:** Each test device includes internal controls to determine whether the VerifyNow IIb/IIIa Test result is valid and causes the instrument to report an error message for an invalid test. The negative internal control test is performed during the first five seconds of the assay, before platelet activation and binding to the fibrinogen-coated beads occurs.

Wet Quality Control

VerifyNow WQC Level 1 and Level 2 are formulated at clinically relevant levels and can be used as part of a laboratory quality control program. The VerifyNow WQC Level 1 and Level 2 acceptable ranges are printed on the test device pouch.

The manufacturer recommends that a Level 2 WQC be run once each time a new lot or a new shipment of VerifyNow IIb/IIIa Test kits is received. VerifyNow WQC (Catalog #85047) is available from the manufacturer. For more information, refer to the WQC package insert.

Call Accumetrics Technical Support for questions regarding WQC materials.

TROUBLESHOOTING

Under certain conditions, the instrument may display an Error or Alarm message. Refer to the VerifyNow System User Manual for a more detailed explanation of these messages. For additional troubleshooting, contact Accumetrics Technical Support at: (telephone) +1 (858) 643-1600; (e-mail) support@accumetrics.com.

CALIBRATION

VerifyNow IIb/IIIa Test devices are calibrated at the factory. This calibration information is contained in the barcode located on each test device, and is automatically read by the instrument.

INTERPRETATION OF RESULTS

The VerifyNow IIb/IIIa Test reports patient results in Platelet Aggregation Units (PAU), which are calculated as a function of the rate and extent of aggregation. The VerifyNow instrument also allows the user to calculate the percent inhibition for a given patient provided that a baseline PAU value was obtained prior to administration of drug. The use of GP IIb/IIIa inhibitors such as abciximab and eptifibatide has been shown to be useful as an antithrombotic strategy.¹ In patients receiving abciximab or eptifibatide and in animal models of coronary thrombosis with abciximab, the ex vivo platelet aggregation response was virtually abolished when GPIIb/IIIa blockade was maintained at steady state >80%.^{2,3} Further study has shown that with all GPIIb/IIIa inhibitors, including abciximab and eptifibatide, platelet function inhibition >95% at 10 minutes after the start of therapy was associated with a significant decrease in the incidence of a major adverse cardiac event.⁴

Reference Ranges

The VerifyNow IIb/IIIa Test was performed in clinical studies at 9 centers on patients undergoing percutaneous coronary intervention. Testing was performed on samples from 173 patients with planned treatment with abciximab (173 baseline and 173 post abciximab) and 114 patients with planned treatment with eptifibatide (114 baseline and 97 post eptifibatide). **Table 2** presents the mean PAU values at Baseline (prior to abciximab or eptifibatide administration), >80% and >95% inhibition.

Abciximab				Eptifibatide		
PAU	Baseline	≥80%	≥95%	Baseline	≥80%	≥95%
n	173	173	130	114	97	50
Mean	224	9	2	199	11	3
SD	54	20	4	40.6	9.4	3
Reference Range	125-330	0-44	0-13	136-288	0-31	0-10

Table 2: Reference Ranges

LIMITATIONS

Delays in testing or difficulty of specimen collection may result in spurious values.

The lyophilized reagent is hygroscopic and can degrade after prolonged exposure to room air. Therefore, the test device should be used shortly after removal from the foil pouch.

When results are not within the expected limits, the possibility of improper sample collection or handling should be investigated. Repeat the test using a new test device and sample. Patients with inherited platelet disorders such as von Willebrand Factor Deficiency, Glanzmann Thrombasthenia and Bernard-Soulier Syndrome have not been studied with the VerifyNow Ilb/IIIa Test. The VerifyNow Ilb/IIIa Test is not intended for use with these types of platelet disorders. Patients with a known history of platelet counts <100 x10⁹/L have not been studied. The VerifyNow Ilb/IIIa Test has not been compared to platelet aggregometry in patients with a 30% decrease in PRP platelet count. VerifyNow Ilb/IIIa Test results should be interpreted in conjunction with other clinical and laboratory data available to the clinician.

Reference ranges were calculated using CLSI guidelines. The following demographic characteristics were considered: gender, ethnicity, current smoker, and previous abciximab use. The mean baseline PAU was consistent across all groups, and no statistically significant bias was demonstrated in any of the demographic variables tested. An analysis of the clinical data and the associated demographic information showed that there was no need to create separate reference ranges based upon demographic variables. As has been shown, there is no overlap between the baseline and inhibition reference ranges, which provides the clinician with an unambiguous interpretation.

PERFORMANCE CHARACTERISTICS

Studies evaluating the VerifyNow Ilb/IIIa Test in patients treated with abciximab compared the VerifyNow Ilb/IIIa Test to a 125 I-abciximab Receptor Blockade Assay (RBA) which quantitates the molecules bound per platelet and platelet aggregometry using a CHRONO-LOG Optical Aggregometer with 20 µM ADP as the agonist to detect the maximal slope over a five minute period. The unit of measure for the VerifyNow Ilb/IIIa Test is PAU and is based upon both the rate and extent of platelet and fibrinogen coated beads binding over a defined period.

Studies evaluating the VerifyNow Ilb/IIIa Test in patients treated with eptifibatide compared the VerifyNow Ilb/IIIa Test to the same platelet aggregometry test described above and a platelet flow cytometric test measured as mean fluorescence units.

The clinical trials were designed to obtain samples at three time points: 1) Baseline, prior to abciximab or eptifibatide administration; 2) During Infusion, within 1 hour following bolus administration and 3) Post Infusion, 24 hours after bolus or at the time of discharge, whichever occurred first, to assess recovery of platelet function. For each patient, a baseline, during infusion and post infusion (24-hour) sample was measured by each test method and the results recorded.

Figures 2 and 3 show the time course of platelet inhibition for the three methods as individual points and mean +/- standard deviation, respectively, and illustrate the agreement of the three methods. (See Figures 2 and 3)


Figure 2: Mean and Standard Deviation of Percent Baseline GPIIb/IIIa or Aggregation (all methods) as a function of Time after Abciximab Infusion


Figure 3: Mean Percent Baseline GPIIb/IIIa Aggregation (all methods) as a function of Time after Eptifibatide Infusion

Method Correlation

Correlation between methods was assessed by regression analysis based on data from patients treated with abciximab and eptifibatide. The data was analyzed to determine the correlation of the VerifyNow IIb/IIIa Test with RBA and platelet aggregometry in abciximab treated patients and with PAC1 Flow Cytometry and platelet aggregometry in eptifibatide treated patients. Data points for the clinical analysis were reported in raw units for each test: PAU for the VerifyNow IIb/IIIa Test; slope for aggregometry; molecules per platelet bound for RBA; and mean fluorescence for PAC1 Flow Cytometry. (See Figures 4, 5, 6 and 7)


Figure 4: Deming Regression VerifyNow IIb/IIIa Test vs Aggregometry


Figure 5: Deming Regression VerifyNow IIb/IIIa Test vs Receptor Blockade


Figure 6: VerifyNow IIb/IIIa Test vs Aggregometry Eptifibatide


Figure 7: VerifyNow IIb/IIIa Test vs PAC1 Flow Cytometry Eptifibatide

Analytical Sensitivity

The VerifyNow IIb/IIIa Test was run (n=6 per patient) on sixteen patients dosed with abciximab and measured at 10 minutes post bolus dose. The analytical sensitivity of the test was determined to be 8 PAU (2 X SD).

Linearity

Figure 8 shows the linearity of the VerifyNow System. The testing consisted of multiple donors ex-vivo dosed at nine eptifibatide concentrations over multiple days with multiple replicates. The linearity data show that the relationship between eptifibatide concentration and % inhibition is not linear.


Figure 8: Linearity of VerifyNow System n=29, ex-vivo dosed, multi-day

Precision

Simple and complex precision for both abciximab and eptifibatide were calculated according to CLSI precision guidelines using three lots of VerifyNow IIb/IIIa Test devices and two instruments for abciximab and two lots of test devices and three instruments for eptifibatide. Whole blood was drawn from volunteer donors over 20 days. The within-run CV shown in Tables 2 and 3 is defined by CLSI as the within-run precision value. (See Tables 3 and 4)

Complex Precision Whole Blood				Within-Run (Simple)		Total	
Description	n (days)	n (results)	Mean	SD	CV %	SD	CV %
Male Age 32	20	78	147.8	9.2	6.2	13.9	9.4
Male Age 28	20	77	185.3	12.2	6.6	16.0	8.6

Table 3: Complex Precision Statistics (Eptifibatide) as Defined by CLSI

Complex Precision Whole Blood				Within-Run (Simple)		Total	
Description	n (days)	n (results)	Mean	SD	CV %	SD	CV %
Female Age 26	20	80	186.3	9.8	5.2	18.7	9.9
Female Age 41	20	80	187.0	9.1	4.9	26.1	13.9
Male Age 29	20	80	244.8	11.1	4.9	25.9	11.5

Table 4: Complex Precision Statistics (Abciximab) as Defined by CLSI

Simple precision was determined for the VerifyNow IIb/IIIa Test wet quality control Level 2. Three lots of test devices were each tested 20 times using a single lot of control. This study was later repeated with a new control lot. The results are shown in Table 4. (See Table 5)

Lot to lot Cartridge Variation				
Device Lot	n	PAU Mean	SD	CV*
1	20	202	12	6.0
2	20	180	15	8.3
3	20	176	17	9.9

Table 5: Simple precision

*The manufacturer's specification for the coefficient of variation is <10%.

Interfering Substances

Clinical study data were evaluated to assess the effect of potentially interfering medications on PAU values. Baseline PAU values were compared for patients with and without the following medications present, and no differences were observed: aspirin, ticlopidine, clopidogrel, heparin, warfarin, acetaminophen, NSAIDs, beta-adrenergic blockers, Ca++ channel blockers, statins and nitrates.

Aspirin and heparin were further studied because they are commonly coadministered with GP IIb/IIIa therapy. Varying concentrations of both drugs were evaluated in in-vitro studies, and no significant effect on PAU values was observed. These data are consistent with data reported by Collier using an earlier prototype version of the VerifyNow IIb/IIIa Test.⁵

The effect of platelet counts, hematocrit values and fibrinogen concentration were studied in two clinical trials. Test performance was not affected by platelet count values between 98,000 and 376,000 platelets per microliter, or hematocrit values between 23 and 45%. No test interference was observed when samples with fibrinogen levels between 104 and 617 mg/dL were tested with the VerifyNow IIb/IIIa Test. Moderate lipemia was found to induce an increase in PAU values by 16%.

Another intravenous GP IIb/IIIa inhibitor, tirofiban (Aggrastat®) is approved for use in patients with unstable angina. Patients who have been treated with tirofiban may be sent to the cardiac catheterization lab for PCI. The VerifyNow IIb/IIIa Test can detect platelet inhibition by this drug; therefore VerifyNow IIb/IIIa Test results for these patients should be interpreted with care. A sample taken prior to abciximab (ReoPro®) or eptifibatide (Integrilin®) administration cannot be used to establish a "baseline" or uninhibited result if a GP IIb/IIIa inhibitor was administered within the past 10 days.

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BIBLIOGRAPHY

1. Collier BS. Blockade of Platelet GPIIb/IIIa Receptors as an Antithrombotic Strategy. *Circulation*. 1995; 92:2373-2379.
2. Mascelli MA, Worley S, Veriabo NJ, Lance ET, Mack S, Schaible T, Weisman JF, Jordan RE. Rapid Assessment of Platelet Function With a Modified Whole-Blood Aggregometer in Percutaneous Transluminal Coronary Angioplasty Patients Receiving Anti-GPIIb/IIIa Therapy. *Circulation*. 1997; 96:3860-3866.
3. Tcheng JE, Talley JD, O'Shea JC, Gilchrist IC, Kleiman NS, Grines CL, Davidson CJ, Lincoff AM, Califf RM, Jennings LK, Kitt MM, Lorenz TJ. Clinical Pharmacology of Higher Dose Eptifibatide in Percutaneous Coronary Intervention (The PRIDE Study). *The American Journal of Cardiology*. 2001;88: 1097-1102.
4. Steinhubl SR, Talley JD, Braden GA, Tcheng JE, Casterella PJ, Moliterno DJ, Navetta FI, Berger PB, Popma JJ, Dangas G, Gallo R, Sane DC, Saucedo JF, Jia G, Lincoff MA, Theroux P, Holmes DR, Teirstein PS, Kereiakes DJ. Point-of-Care Measured Platelet Inhibition Correlates With a Reduced Risk of an Adverse Cardiac Event After Percutaneous Coronary Intervention: Results of the GOLD (AU-Assessing Ultegra) Multicenter Study. *Circulation*. 2001;103:2572-2578.
5. Collier BS, Lang D, Scudder LE. Rapid and Simple Platelet Function Assay to Assess Glycoprotein IIb/IIIa Receptor Blockade. *Circulation*. 1997; 95:860- 867.
6. Kini AS, Richard M, Suleman J, Perez N, Lee P, Fisher EA, Kamran M, Marmur JD, Sharma SK. Effectiveness of Tirofiban, Eptifibatide, and Abciximab in Minimizing Myocardial Necrosis During Percutaneous Coronary Intervention (TEAM Pilot Study). *The American Journal of Cardiology*. 2002;90:526-529.
7. Phillips DR, Teng W, Arfsten A, Nannizzi-Alaimo L, White M, Longhurst C, Shattil S, Randolph A, Jakubowski J, Jennings L, Scarborough R. Effects of Ca²⁺ on GPIIb/IIIa Interactions with Integrilin. *Circulation*. 1997;96: 1488- 1494.

EXPLANATION OF SYMBOLS / SYMBOLES / LEGENDA DEI SIMBOLI /
ERLÄUTERUNG DER SYMBOLE / EXPLICACIÓN DE LOS SÍMBOLOS

	Catalog Number Numéro de catalogue Numero di catalogo Katalognummer Número de catálogo
	Origin-Human Origine : humaine Origine: umana Ursprung: human Origen humano
	In vitro Diagnostics Device Pour diagnostics in vitro Dispositivo diagnostico in vitro In-vitro-Diagnostikum Dispositivo diagnóstico in vitro
	Manufacturer Fabricant Produttore Hersteller Fabricante
	Temperature Limits Limites des températures Limiti di temperatura Temperaturbereich Límites de temperatura
	Use By Utiliser avant Da utilizzare entro Verfallsdatum Fecha de caducidad
	Authorized Representative in the European Community Représentant agréé au sein de la Communauté européenne Rappresentante autorizzato per l'Unione Europea Autorisierte Vertretung in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea
	Lot Number Numéro de lot Numero di lotto Chargennummer Número de lote
	Contains sufficient for < n > tests Possibilité de réaliser < n > tests Contiene materiale sufficiente per < n > analisi Ausreichend für < n > Tests Contiene cantidad suficiente para < n > pruebas
	Caution, consult accompanying documents. Attention, consulter les documents joints Attenzione, consultare la documentazione allegata. Achtung, Begleitdokumente einsehen. Precaución, consultar la documentación adjunta.
	Consult Instructions for Use Consulter le mode d'emploi Consultare le Istruzioni per l'uso Gebrauchsanweisung einsehen Consulte el modo de empleo

Accumetrics, Inc.

3985 Sorrento Valley Blvd.

San Diego, California 92121 U.S.A.

Telephone: +1 858 643 1600

Email address: support@accumetrics.com

Website: www.accumetrics.com

MDSS GmbH

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30175 Hannover, Germany